




Tips for ENT

SS2

High-effectiveness osteotomy.




Model: Bone



SS1L

Left angle 90°, osteotomy.




Model: Bone



SS1R

Right angle 90°, osteotomy.




Model: Bone



SS5

High-efficiency osteoplasty.




Model: Bone



SS8

Multi-faceted distribution of blade teeth, suitable for various osteotomies.




Model: Bone



SS11L

Left angle 90°, circumferential teeth, for osteotomy in special bone positions.




Model: Bone



SS11R

Right angle 90°, circumferential teeth, for osteotomy in special bone positions.




Model: Bone





SS13

For osteotomy.




Model: Bone



SS19L

Left curved, suitable for oblique cutting in osteotomy.




Model: Bone



SS19R

Right curved, suitable for oblique cutting in osteotomy.




Model: Bone



SS21

For osteotomy: fan-shaped blade design for easy multi angle cutting.




Model: Bone



SP7

Diamond-coated ($40\mu\text{m}$), Used for safe bone perforation. Please note that it can only be used in Perio mode.




Model: Perio



SP8

Coarse mesh-shaped concave-convex surfaces on both sides of the head, for osteoplasty and bone trimming. Please note that it can only be used in Perio mode.




Model: Perio



SP9

Finely mesh-shaped concave-convex surfaces on both sides of the head, for osteoplasty and bone trimming. Please note that it can only be used in Perio mode.




Model: Perio



SS35

The head is in the shape of a spatula, used for osteoplasty and bone repair.




Model: Bone





SS36

Osteotomy.




Model: Bone





SS37

Osteotomy.




Model: Bone





SS38

Used for transverse bone cutting in osteotomy.




Model: Bone





SP18

The head plane is rough and concave convex, used for osteoplasty and bone repair surgery.




Model: Perio



SP19

The head plane is rough and concave convex, used for osteoplasty and bone repair surgery.



Model: Perio

